

Fit Guide-Bantam Extra Small and Small-Read all instructions thoroughly before beginning any adjustments.

Foot Plates

Measure the distance from the bottom of the foot to the bottom of the leg behind the knee. (A)

Before the child is in the stander, loosen the clamps (do not remove) on the foot plates and adjust to the height that was measured for the foot plate to seat height. (A) Tighten securely.

⚠ **Caution:** Do not adjust the seat depth while the child is in the unit.

⚠ **Caution:** Ensure the colors on the seat tube and following arm match.

Measure the distance from the back of knee to the back of the buttocks. (B)

Remove the Shadow Tray and bring the seat up to a slight angle. (B)

Use the following steps to adjust the seat depth. Use the measurement B from the front of the seat to the seat back. (B)

1. Loosen the knobs on the seat tube and the following arm. Do not unscrew knobs out all the way when adjusting the seat. (C)

The **Bantam Extra Small** following arm (back angle) adjustments are located above the locking knob. The back angle adjustments go from green (largest) to blue (smallest). (C)

The **Bantam Small** following arm (back angle) adjustments are located below the locking knob. The back angle adjustments go from blue (smallest) to green (largest). (C)

2. **Adjusting to a shorter setting:** Adjust the following arm (back angle) first by depressing the spring button and moving no more than two holes at a time towards the blue setting. (F)

3. Adjust the seat tube (seat depth) one hole at a time. Continue this sequence until you achieve the desired setting and the colors on the seat tube and following arm match, towards the blue setting. (D)



4. **Adjusting to a longer setting:** Adjust the seat tube (seat depth) first by depressing the spring button and moving only one hole at a time. (D)

5. Adjust the following arm (back angle) no more than two holes at a time. Continue this sequence until you achieve the desired setting and the colors on the seat tube and following arm match. (E)

Note: For easier adjustment you may need to lift the seat slightly as the seat depth adjustment is made.

Tighten the seat depth and following arm knobs securely when you have achieved the desired position.

Inspect your standing frame periodically to ensure it is in safe operating condition



Fit Guide-Bantam Medium-

Read [all instructions](#) thoroughly before beginning any adjustments.

Measure

Before transferring into the Bantam medium, measure for approximate seat depth. Measure from the back of the knee to the back (A). Measure the seat to foot plate settings by measuring from the bottom of the leg behind the knee, to the bottom of the heel (B).

If the supine option was ordered, make sure the controller is in the “Sit-to-Stand” position prior to adjusting seat depth (C). The seat should be lowered to a seated position for easiest adjustment.

Seat Depth

Use the measurement (A) to adjust the seat depth on the stander, measure from the front of the seat to the back.

To adjust the seat depth, loosen the knob under the seat (D). With one hand supporting the back (if equipped) depress the spring button on the outer left side of the seat tube underneath the seat (E). Using the back as a lever, move the seat depth to the desired setting, noting the color. Tighten knob securely (D).

Be sure to have one hand supporting the back assembly (F). Loosen the knob on the following arm (G). With one hand supporting the back assembly flip the gray plunger lever in the center of the knob fully open (H). Then continue to loosen the knob until resistance is felt. DO NOT continue to turn the knob counterclockwise. Forcing the knob past its stop will damage the mechanism.

⚠ Caution: Make sure to have one hand supporting the back assembly, flipping the gray plunger lever will cause the back assembly to move freely.

Using the back as a lever, adjust the following arm so that the color of the seat depth decal matches the color on the following arm in the indicator window. Release to close the gray plunger lever to re-engage the locking plunger in the clamp knob. You may need to partially tighten the knob and slightly move the back assembly to allow the locking pin to engage the “home” position. Tighten the knob securely.

Back Angle



After achieving the desired seat depth, rotate the back angle adjustment knob (I) so the seat/back is in the desired angle.

Foot Plates

Use the (B) measurement to adjust the foot plate height. Loosen the ratchet handles (J) on the foot plates to adjust height and plantar/dorsi. Tighten ratchet handles after positioning. The knees should be slightly lower than the hip with the legs resting on the seat.

Knee Pads

Position the user’s hips as close as possible to stander’s hip pivot (K). Position the knee pads over the knees, not leaving more than a finger space between the back of the knee and the seat (L). Tighten knees securely (M).



Inspect your standing frame periodically to ensure it is in safe operating condition

Fit Guide-Evolv Medium, Large and XT - Glider Medium and Large

Read all instructions thoroughly before beginning any adjustments.

Measure

1. Before transferring into the Evolv, measure for approximate (A) seat depth and seat to (B) foot plate settings.

⚠ **Caution:** These adjustments are only to be made while the unit is unoccupied.

Seat Depth

2. Use the above measurement (A) to adjust the seat depth on the stander, measure from the front of the seat to the back. To adjust the seat depth, loosen the knob under the seat (C). Depress the snap button on the adjustment tube while moving the back (D). Adjust the seat depth until the proper measurement has been reached. Tighten the knob (C) after adjusting.

⚠ **Caution:** The following adjustments are only to be made in the seated position.

Back Angle

3. After achieving the desired seat depth, rotate the back angle (E) adjustment knob so the (F) seat/back angle indicator is in the green area prior to standing.

Note: Repeat steps 2 and 3 for each color adjustment.

Foot Plates-Evolv

4. Use the (B) measurement to adjust the foot plate height. Loosen the (G) ratchet handles on the foot plates to adjust height and plantar/dorsi. Tighten ratchet handles after positioning. The knees should be slightly lower than the hip (H).

Foot Plates-Glider

5. Use the (B) measurement to adjust the foot plate height. Loosen the (G) knobs on the foot plates to adjust height. Tighten the knobs after positioning. The knees should be slightly lower than the hip (H).

Positioning

5. Position the user's hips as close as possible to (I) Evolv's hip pivot. To determine correct positioning on the Evolv, only two fingers should fit between the front of the seat and the back of the knee (J).

Knee Pad-Evolv


6. Adjust the (K) knee pad so it fits comfortably up to the knees.

Knee Pad-Glider


6. Slide into knee pads (J) and secure the calf straps behind the knees.




(A)



(B)




(C)



(D)



(E)



(F)

(G) Evolv

(G) Glider



(H)



(I)



(J) Evolv



(K) Evolv



(K) Glider




Inspect your standing frame periodically to ensure it is in safe operating condition

Fit Guide-StrapStand-*Read all instructions thoroughly before beginning any adjustments.*

1. To use the StrapStand, position the wheelchair as close as possible to the knee pad. Casters on the StrapStand should be locked.
2. Position feet in the foot plates. Adjust knee pad height so knees are centered in the knee pad. With the wheelchair locked, slide forward on the seat so that the knees make contact with the knee pad.
3. Hook the strap to one of the lifting arms and bring the other side around the body. Then hook the strap to the other lifting arm.
4. After attaching the strap to the lifting arms, adjust the hook position. The strap can be adjusted to accommodate different size individuals. This will vary depending on the strap size and the girth of the individual.
5. Adjust the placement of the strap under the buttocks so that the top edge of the strap is slightly under the belt line. Slide it into place as shown.
6. Partially lift the user by pumping the actuator handle towards the user to tension the strap. At this time and if necessary, reposition the strap. Continue until the user is in the standing position. The user's wheelchair should remain in position and available so the user can return to the seated position.

⚠ **Caution:** Never adjust the strap while in the standing position.



Fit Guide-Zing Size 1 and Zing Size 2

Read all instructions thoroughly before beginning any adjustments.

Centering the Hip (A):

The user's hip should be centered in line with the pivot point of the stander ensuring proper hip alignment when using abduction, adjust support components up or down based on hip position.

Foot Plate Height (B):

Measure the distance from the hip of the user to the bottom of the foot.

Knee Pad Height (C):

Measure the distance from the hip of the user to center of the knee.

