



Taking you where you want to be

Discover the  
easiest way to  
add power  
assist to your  
manual  
wheelchair!



Introducing AMP – the innovative device that brings power assist to your manual wheelchair – all without attaching permanent hardware, motors or wheels!

- **Rolls into place:** Simply roll into the AMP and it automatically clicks into position on your chair – no modifications required.
- **Starts with just a push:** Just give your wheelchair's handrims an easy push to get the AMP moving.
- **Navigates naturally:** Use your existing manual handrims to start, stop and change directions – just as you do now but with less effort.
- **Responds to your commands:** Easily power up or change mode settings with the AMP's intuitive Bluetooth remote.
- **Protects your shoulders, wrists & hands:** The AMP gives you a little extra help when you need it.
- **Undock when you're done:** When you don't need the AMP, simply roll off it and use your wheelchair like always with no extra weight to push around.



Get AMP'd!

Proudly distributed in Australia by







*"The AMP is my superpower – it's the power assist I've been waiting for! I feel like I'm living life again and not missing out on the little things like enjoying outings and walking the dogs!" - Tom Mullen*



## AMP Power Assist with Handrim Input



Easy roll-on, roll-off



Bluetooth remote



Optional joystick input

## AMP Handrim Specifications

- Available for: Most rigid or folding wheelchairs with 22", 24", 25" or 26" wheels. Some restrictions apply.
- Wheelchair width compatibility: 14"-20"
- Bluetooth remote
- Optional joystick
- Power module weight: Less than 42 lbs. with battery.
- Battery weight: 6.5 lbs.
- Weight limit: 300 lbs.
- Unique Life P04 Lithium Iron Phosphate battery – safest battery technology in lithium family.
- Battery charging: Standard charger: 6.5 hours, Super charger 2-2.5 hours.
- Range: Up to 11 miles on single charge.
- Speed: Up to 4 mph.
- HCPC Code: E0986

# Get AMP'd!



P: 1300 561 334  
E: [sales@paragonmobility.com.au](mailto:sales@paragonmobility.com.au)  
W: [paragonmobility.com.au](http://paragonmobility.com.au)

MM-AMPF 2.25