

Indoor Bases for Stroller Seat Units

Height-adjustable indoor bases are available for the BINGO Mini, BINGO Evo, and BINGO Spex, where the removable seat units can be removed from the stroller base and attached to the indoor base with ease. Height-adjust features are quick and easy to change allowing the height of the seat to be lowered to the floor for peer interaction, raised to the dinner table to make mealtimes enjoyable for the whole family, or raised even further to breakfast bar height for easier feeding!

Cobra Base

Height-adjustable indoor base to suit the **BINGO Mini** (see page 14), **BINGO Evo** (see page 16), and **BINGO Spex** (see page 18) seat units.

German designed precision technology with a focus on quality and durability

45° of tilt

Simple and quick
Hand Operated Gas Assist height adjustment



Compact footprint for easy manoeuvrability and navigation around table legs

Lockable castors to stop the base from rolling

Size	Code	User Age	Max Load	Approx Product Weight	Tilt	Seat Height	Total Length	Total Width
1	2880-8711-000	1-5 years	50 kg	11.3 kg	45°	27-61 cm	71 cm	61 cm
2	2880-8721-000	3-15 years	50 kg	11.5 kg	45°	27-61 cm	71 cm	61 cm