

# BINGO Evo Stroller

Larger supportive seating system with many features, including a wide range of accessories and options to meet the needs of different children.

## Modular seating & positioning stroller

Multi-directional suspension **protects the child** by absorbing ground bumps, and is effective off road or on rough terrain

**Easily adjust to Rear Facing configuration** for improved close interaction with caregiver or monitoring of safety

**Stepless seat depth adjustment**

**Note:** specialty cover not included, must be added as an accessory.



**Ergonomic,** height-adjustable push handle

Sun canopy option available to allow parents and caregivers to still get out and about, knowing the child is protected from the sun

**Adjustable seat unit**  
Seat depth and width, back support height and lower leg length

### Evo Package Includes:

Removable seat, matte black mobility base with swivel front black mag wheels and solid tyres, tie-downs, hip guides, basket, adjustable suspension and 3 sided foot support.



### Colour Options:



Turquoise Soft (-349)

Berry Soft (-347)

Black (-346)

### ONE SEAT. TWO BASE OPTIONS

Indoor Cobra base option available. Height adjustable while the child is seated.



### ADAPTIVE SUSPENSION SYSTEM



Fourfold adjustable frame suspension system which can be easily adjusted without tools.

### TILT

The special construction of the tilt mechanism allows a 45° tilt of the seat unit without compromising the overall stability.



### COMPACT FOLDED SIZE



Easier to store in a car boot! The innovative folding mechanism can be released easily on the push bar.

### BACK SUPPORT RECLINE

The operating elements of the tilt & back support recline are positioned centrally and allow for easy and secure handling.



### NEW TUBIE FUN COVERS

Personalise and protect your Bingo Evo Stroller with our brand-new, waterproof and washable covers from Tubie Fun!



Size	Code	Seat Load Capacity	Overall Length	Overall Width	Seat Depth	Seat Width	Back Support Recline	Shoulder Height	Foot Support Height	Overall Weight	Indoor Base Seat Height	Indoor Base Footprint
1	2881-0008-000	50 kg	92 cm	63 cm	16-34 cm	18-30 cm	90-160°	28-46 cm	16-37 cm	15 kg	27-67.5*	61cm (W) x 71cm (L)
2	2882-0008-000			68 cm	22-40 cm	23-35 cm	90-160°	34-55 cm	16-40 cm	16.1 kg		
3	2883-0008-000			68 cm	27-45 cm	23-35 cm	90-135°	42-63 cm	16-40 cm	18 kg		

\* Depending on size of castors

# Accessories

**Occipital Head Support**  
Merino/Spacer Fabric  
2010-4110-099  
**Wipe Down**  
2010-4110-209



**Std Lateral Head Support,**  
Merino/Spacer Fabric  
Small: 2010-4210-099  
Large: 2010-4211-099



**Standard Lateral Head Support, Wipe Down**  
Small: 2010-4210-209  
Large: 2010-4211-209



**Ext Lateral Head Support,**  
Merino/Spacer Fabric  
Small: 2010-4310-099  
Large: 2010-4311-099



**Extended Lateral Head Support, Wipe Down**  
Small: 2010-4310-209  
Large: 2010-4311-209



**Std Lateral Trunk Support**  
1: 2880-3010-346  
2: 2880-3020-346  
3: 2880-3030-346



**Spex Lateral Trunk Supports**  
Narrow Straight: 2880-3008-100  
Broad: 2880-3008-310  
Teardrop: 2880-3008-500



**Additional Padding**  
Sacral Support: 2020-2510-060  
Thigh Wedge: 2006-0505-000



**Medial Knee Support**  
1: 2880-1651-346  
2: 2880-1652-346  
3: 2880-1653-346



**Anterior Chest Support**  
1: 2880-6751-346  
2: 2880-6752-346  
3: 2880-6753-346



**Wonderseat Padded Centrepoint Harness**  
Size 0: 1409-5660-017  
Size 1: 1409-6661-017



**5-Point Anterior Chest Support**  
Bingo: 2880-6511-346  
Spex: 2006-6551-000



**Hip Belt**  
Spex 2-Point: 1405-6220-027  
Wonderseat 4-Point: 1406-5420-027



**Groin Strap**  
1: 2880-6361-346  
2-3: 2880-6362-346



**Foot Support Spacer 4cm thick**  
1: 2880-5321-000  
2: 2880-5322-000



**Ankle Foot Fast**  
1: 2880-6101-000  
2: 2880-6102-000



**Therapy Tray**  
1: 2880-7111-000  
2-3: 2880-7112-000  
Infill: 2880-7180-000



**Grab Rail**  
1: 2880-1651-346  
2-3: 2880-1652-346



**IV Pole:**  
2880-8275-000  
**Extra Bracket**  
2880-8275-050



**Oxygen Bottle Holder:**  
2880-8253-000



**Chassis Mount Tray**  
1: 2880-8058-000  
2: 2880-8059-000



**Separate Seat Unit**  
1: 2881-0011-000  
2: 2882-0012-000  
3: 2883-0013-000



**Sun Canopy & Rain Shield:**  
2880-7500-346



**Cobra Indoor Base**  
1: 2880-8711-000  
2: 2880-8721-000



**Head Support Pads**  
1 2880-4471-346  
2 2880-4472-346



**Swivel Lock:**  
2880-8440-000



**Cover Set - Llamas**  
Size 1: 2880-9561-121  
Size 2: 2880-9562-121



**5 Point Harness Cover:**  
2880-6511-121

**Cover Set - Planets**  
Size 1: 2880-9561-036  
Size 2: 2880-9562-036



**5 Point Harness Cover:**  
2880-6511-036