

Bingo Evo Mini Stroller

Medifab's newest early intervention special needs stroller with a wide range of accessories and options to meet many requirements.

Specialised Early Intervention Tilt-in-space Stroller



Adjustable seat unit
Seat depth and width, backrest height and lower leg length

Ergonomic,
height-adjustable push handle

Individual upholstery
Removable seat, back, head and laterals cushions allow easy customisation

Multi-directional suspension **protects the child** by absorbing ground bumps, and is effective off road or on rough terrain

Stepless seat depth adjustment

Sun canopy option available to allow parents and caregivers to still get out and about!

Mini Package Includes:

Removable seat, matte black mobility base with swivel front black mag wheels and solid tyres, tie-downs, hip guides, basket, adjustable suspension, receiver for supports & footplate.



Colour Options:



Turquoise (-339) Pink (-340) Red (-345)

Size	Code	Seat Load Capacity	Push Handle Height	Overall Width	Seat Depth	Seat Width	Backrest Recline	Shoulder Height	Footplate Height	Overall Weight	Indoor Base Seat Height	Indoor Base Footprint
0	2880-0008-000	25 kg	75-110 cm	63 cm	12-27 cm	13-25 cm	80-180°	25.5-39.5 cm	16-34 cm	17.6 kg	27-67.5*	61cm (W) x 71cm (L)

* Depending on size of castors

ONE SEAT. TWO BASE OPTIONS

Indoor Cobra base option available. Height adjustable while the user is seated.



ADAPTIVE SUSPENSION SYSTEM



Fourfold adjustable frame suspension system which can be easily adjusted without tools.

TILT-IN-SPACE

The special construction of the tilt mechanism allows a 45° tilt of the seat unit without compromising the overall stability.



COMPACT FOLDED SIZE



Easier to store in a car boot! The innovative folding mechanism can be released easily on the push bar.

BACK SUPPORT RECLINE

The operating elements of the tilt & backrest recline are positioned centrally and allow an easy and secure handling.



Accessories

<p>Head Support Pads 0: 2880-4470-323</p> 	<p>Std Lateral Support 0: 2880-3009-323 Pressure-relief upgrade 2880-3846-400</p> 	<p>Occipital Head Support Spacer: 2010-4101-099 Wipe Down: 2010-4101-209</p> 	<p>Std Lateral Head Support Spacer: 2010-4201-099 Wipe Down: 2010-4201-209</p> 	<p>Ext Lateral Head Support Spacer: 2010-4301-099 Wipe Down: 2010-4301-209</p> 
<p>Upholstery pressure-relief foam upgrades Seat: 2880-0391-000 Back: 2880-2392-000</p> 	<p>Seat options Sacral Support: 2020-2510-060</p> 	<p>Seat options Thigh Wedge: 2006-0505-000</p> 	<p>Seat options Trochanter Wedges: 2006-0660-000</p> 	<p>Abductor Pommel 1: 2880-1651-323</p> 
<p>Footplate Spacer 2880-5321-000</p> 	<p>Neoprene Footplate Cover 2880-5305-000</p> 	<p>Ankle Foot Fast 2880-6101-000</p> 	<p>Therapy Tray Can't be used with grab rail 1: 2880-7111-000 Infill: 2880-7180-000</p> 	<p>Grab Rail Can't be used with tray 2880-7250-323</p> 
<p>Oxygen Bottle Holder: 2880-8253-000</p> 	<p>Chassis Mount Tray 2880-8058-000</p> 	<p>Sun Canopy: With mounts & rain cover 2880-7501-323</p> 	<p>Accessory Bag Handle-bar strap system 2880-8012-099</p> 	<p>Hand Strap: For AS/NZS 3696.19 2005-6930-000</p> 
<p>Swivel Lock: 2880-8440-000</p> 	<p>Cobra Indoor Base 0: 2880-8720-000</p> 	<p>Spex 5-Point Harness: 2006-6551-000</p> 	<p>IV Pole: 2880-8275-000</p> 	