

Carrot 3000 Car Seat



Various extension accessories and multi-adjustable features allow for growth and adaptation as a child's needs change over time. The Carrot 3000 caters for safety and the importance of keeping the child in a supported posture for travel.

Our most comprehensive special needs car seat

Basic package includes

- Carrot seat unit
- 5-point anterior chest support
- Headrest wings
- Tee-pad
- Free recline back support

Ergonomically designed **head support** holds the head and protects from side impacts and centrifugal forces

Tether options include additional security tether for the top of back support or ISOfix latch connectors. These may be required by your local authorities

2 positions of **seat depth adjustment** for a better fit



Shoulder protector wings: Extendable wings hold and protect the body from side impacts and centrifugal forces

Large **chest support** holds the body securely

Recline function with angle adjustable seat and back support



Huge growth for pre-schoolers to older children

When requesting quote, please specify the user's shoulder height and seat depth.



Size	Code	Age Range	User Weight Range	Seat Depth	Seat Width	Shoulder Height	Foot Support Height	Seat Weight
Small	4116-0010-259	3-15 years	15-36 kg	25-40 cm	31cm	32-52 cm 45-67 cm * (*when shoulder protector wings are added)	Small = 25-30 cm Medium = 30-40 cm Large = 40-50 cm	8.7 kg

The Carrot 3000 conforms to EU (ECE R44/04) standards.

Large chest pad/head support model available (4116-0000-259)

If user is too large for the Carrot 3000 Car Seat, the Carrot XL Seat (4116-0020-259) is available for larger users. See page 48

Air Travel Friendly
The Carrot 3000 Car Seat has been approved for use in aircraft of some major airlines including the Federal Aviation Administration in the USA. This is NOT a guarantee that your chosen airline will allow the Carrot 3000 Car Seat to be used in their aircraft, and you must contact the airline prior to flying.

Accessories

Swivel Base Attachment*
 4116-0180-104
Swivel Base Protection Plate
 4116-0190-000



* Requires ISOfix

Foot Support
 Sml lower leg (D): 25-30 cm
 4116-5241-250
 Med lower leg (D): 30-40 cm
 4116-5242-250
 Lrg lower leg (D): 40-50 cm
 4116-5243-250



Swivel Base Footplate Adaptor
 4116-5120-000



Seat Belt Holder
 4116-6790-017



Removable Incontinence Cover
 4116-0475-209



**Rodded 50mm
 Seat Depth Pad**
 4116-0281-259



**100mm Rodded Waterfall
 Seat Base Extension**
 4116-0283-259



**50mm Seat Depth
 Growth Pad**
 4116-0285-259



**100mm Seat Depth Growth
 Pad (only with Waterfall)**
 4116-0289-259



Abduction Pad
 4116-1510-000



**50mm Back Support Height
 Extension Infill Pad**
 4116-2481-259



**100mm Back Support Height
 Extension Infill Pad**
 4116-2482-259



Shoulder Protector Wings
 4116-2753-259



**Anti-Escape Anterior
 Chest Support**
 4116-6026-259



Shoulder Strap Retainer
 Small: 4116-6721-000
 Medium: 4116-6722-000
 Large: 4116-6723-000



**Pen Release Anterior
 Chest Support**
 4116-6025-259



**Cross Strap Anterior
 Chest Support**
 Anti-Escape 4604-3004-026
 Pen Release 4604-3004-025
 Side Release 4604-3112-019



Tether Strap
 4116-6953-017



ISOfix Lower Straps
 4116-6910-000



Hinge Lock Belt
 (for easy carrying)
 4116-8170-000



Tray / Forearm Support
 4116-5534-259



Tilt Wedge
 4116-0291-259



**Firm Contoured Lateral
 Trunk Support Wedges**
 4116-3250-090



Air Mesh Cushion Pad
 4116-8880-099

