



EZPZ

Mobile Shower Commode Chairs



EZPZ-AP
Attendant Propel



EZPZ-SP
Self Propel



EZPZ mobile shower commode chairs offer exceptional comfort and function at an affordable price ... but with standard features that are normally found on more expensive chairs!




EZPZ mobile shower commode chairs are exceptionally durable and easy to clean ... ideal for infection-control environments.

FEATURES



Corrosion-resistant aluminum frame with tool-less height adjustment



Height-adjustable without tools (EZPZ-AP) – for toilet bowl clearance and matching transfer heights



Ergonomically contoured seat for best-in-class comfort



Swing-away, removable foot supports – facilitates transfers



Soft and flexible back support – comfortable and easy to clean



Flip-up arm supports – aids in transfers



Dual-locking casters (swivel and roll) – stabilizes chair for toileting, bathing and transfers




Commode pan with lid – includes anti-splash baffle and high-gloss surface to ease cleaning

OPTIONS



Soft seat overlay – molded / contoured foam for added pressure distribution, comfort and positioning



Pelvic belt – side-release buckle (**standard webbing or infection control material**)



Chest belt – 1-piece with side-release buckle (**standard webbing or infection control material**)



Snap-on Infection Control Calf Strap – dries quickly



Arm support locks – provide more security during transfers



Swing-away laterals – provide lateral trunk support



Head support - molded pad with fixture, adjusts in height and fore/aft



Quick-release axles (EZPZ-SP) – makes transporting the chair easier

SPECS

SPECIFICATIONS	EZPZ-AP	EZPZ-SP
Weight Capacity	150 kg (330 lbs)	150 kg (330 lbs)
Seat Width	50 cm (19.5")	50 cm (19.5")
Seat Depth (front of back post to front of seat)	43 cm (17")	43 cm (17")
Back Width	47 cm (18.5")	47 cm (18.5")
Back Height (from top of seat to top of upholstery)	43 cm (17")	43 cm (17")
Toilet Height Clearance:	39-49 cm (15.3"-19.25")	39-49 cm (15.3"-19.25")
Toilet Width Clearance:	43 cm (17")	43 cm (17")
Arm Support Height Above Seat*	23 cm (9")	23 cm (9")
Seat-to-Floor Height* – 2.5 cm (1") increments	48-58 cm (19"-23")	48-58 cm (19"-23")
Overall Width	58 cm (23")	71 cm
Product Weight	16 kg (35 lbs)	18 kg (39.5 lbs)
Frame Material	Aluminum	Aluminum

* without soft seat overlay