

Jenx Multistander

Versatile standing frame, suitable for children from around 9 months right up to 13 years across 2 sizes which provides highly supportive standing. Ideal for growing children in multi-use environments.

Supine, Prone & Upright Standing Frame

De-rotational thoracic & pelvic straps fasten in the centre and adjust independently from either side, enabling a mid-line position to be achieved and maintained. Pads easily remove for side transfers in Supine



Independently adjustable thoracic and pelvic boards cater for the child's individual proportions while giving scope to avoid pressure on delicate areas such as feeding buttons

Deliberately open design ensures you can always see precisely how the user is positioned each and every time

PU boards and support pads are flexible, easy clean and hygienic



Always supportive!

The Multistander is multi-positional meaning it has the flexibility to grow and change with the child easily, supporting them as their skills develop. Prone or Supine configurations mean the correct level of support is always achievable.



Huge growth range. Size 1 shown at smallest height; Size 2 at largest height. (Size 2 provides option for power or manual adjustment).



Angle can be adjusted with the child in the frame and is easily set using the angle gauge.



Set to exactly the right angle with continuous adjustment. Horizontal position assists transfers.

ACCESSORIES



Choose from Lime -054 or Blue -066 colour ways in most accessories

Multigrip Head Support
1 Link: 5356-4601-000
2 Link: 5356-4602-000
3 Link: 5356-4603-000



Moulded Head Support Oval Head Support
5230-4101-XXX
5356-4611-XXX



Slim Thoracic Lateral Trunk Supports
5356-3241-XXX
Hi-Rise Thoracic Lateral Trunk Supports
5230-3251-XXX



Multi-Use Tray
Size 1: 5230-7111-000
Size 2: 5230-7112-000
Elbow Blocks: 5356-7351



Multigrip Body Support Size 2:
5230-2721-XXX



Hip Supports
5356-3221-XXX



Sternal Prompt (Prone)
5230-2910-000
Needs cover 5230-8801



Knee Supports
Size 1: 5356-1751-XXX
Size 2: 5356-1752-XXX



Waistcoat Anterior Chest Support
Size 2: 5230-6712-000
Size 3: 5230-6713-000
Size 4: 5230-6714-000



Sandals
Size 1: 5356-6081-000
Size 2: 5356-6082-000
Size 3: 5356-6083-000



Size	Code	Age	Max User Height	Angle Range	Frame Width	Chest Pad Width	Hip Pad Width	Chest Pad Height	Height of Stander	Foot Support
1	5230-0001	9 mths-6 yrs	116 cm	10° - 90°	50-90 cm	9 cm	66 x 50 x 60 cm	Prone 49-86 cm Supine 57-86 cm	67-94 cm	-15°
2	5230-0002	3-13yrs	175 cm	0° - 80°	67-113 cm	10 cm	86 x 80 x 74 cm	Prone 70-117 cm Supine 77-117 cm	83-128 cm	+15°