

Jenx Multistander with Abduction

The simple abduction accessory for the Multistander size 1 offers children the opportunity to stand in abduction.

Abduction standing for children as young as 9 months

0-30° hip abduction can be achieved in each leg in 5° incremental stages



Abducted standing close to the floor, facilitating play and interaction with peers at eye level

Simple, tool free adjustment using a strong locking mechanism to keep the hip position secure

Offers smaller children (as young as 9 months) the opportunity to continue standing at peer height in either a neutral or abducted hip position

ACCESSORIES

All Jenx Multistander accessories seen on the previous page, are compatible with the Jenx Multistander Abduction accessory.



Retrofit to your Multistander

The Multistander Abduction accessory is available for purchase as a standalone accessory to retrofit to your existing Multistander Size 1.



Code	Angle Range	Age Range	Max User Height	Frame Width	Frame Length	Chest & Hip Pad Width	Knee Dlock Width	Height of Stander	Leg Abduction	Max User Weight	Foot Support Angle Dorsiflexion	Foot Support Angle Plantarflexion
5230-5000-000	0-90°	9 mths - 6 yrs	116cm	56.5cm	71cm	14-27.5cm	12-40cm	61-94cm	30°	35kg	15°	15°

This includes the Abduction Accessory only without the Multistander frame.