



Jumbo prone & upright stander

REHATEC®

The Jumbo upright stander is a standing frame for children's first successful standing training. For children with hip-knee-ankle-foot orthosis, a separate model is available to prevent children from falling. The orthosis is securely attached to the stander.




Jumbo for hip-knee-ankle-foot
orthosis
Frame colour: yellow green



Multi tilting base



Metal loop for
hip-knee-ankle-foot orthosis



Tray with elephant design



Homepage



Product page
Jumbo



GMFCS
Level
III - IV



1 Size



0.65 m -
1.10 m



Max.
40 kg



2 cover
colours



22 frame
colours



Infants



Children

