

# Kidsflex Car Seat

Excellent posture support options and suitable for installation with most car types. Available with a range of options and adjustments, allowing complete comfort and postural management for the child.

## Fully adjustable specialist seating



Secure seat using the **existing 3 point safety belt** fitted in the rear seat.

Additional belt guide can be added for other side of Car Seat  
4170-6960-000

**5 point harness** comes standard with abdominal padding

**Backrest with recline function**  
Increased tilt angle from 10°-20° can be added  
4170-0250-010

**Slimline seat cushions** allow more space and headroom inside the car

### Standard package includes

- Super-slim swivel base (left turning\*)
- Tilt angle adjustment (0-10°)
- Clamp bow

\* Right turning available as special order



Easily adjustable in seat depth.



Recline backrest. Easily adjustable in harness height.



Ample growth opportunities.

## ACCESSORIES

### ISOFix Adapter

- 1: 4170-6911-100  
2/3: 4170-6912-100



### Hip/Thigh Supports

- LH: 4170-1111-000  
RH: 4170-1112-000



### Footrest

- 1: 4170-5241-000  
2/3: 4170-5242-000



### Shoulder Supports

- 4170-2750-000



### Lateral Head Support

- 1: 4170-4201-000  
2: 4170-4202-000  
3: 4170-4203-000



### Extended Head Support

- 1: 4170-4301-000  
2: 4170-4302-000  
3: 4170-4303-000



### Small Lateral Supports

- Fixed: 4170-3301-000  
Swing-away: 4170-3401-000



### Large Lateral Supports

- Fixed: 4170-3302-000  
Swing-away: 4170-3402-000



### Abductor Pommel

- 1: 4170-1511-000  
2: 4170-1512-000  
3: 4170-1513-000



### Support Tray Table

- 1: 4170-5531-000  
2: 4170-5532-000  
3: 4170-5533-000



Size	Code	Max Weight	Seat Width	Seat Depth	Shoulder Height	Lateral Height	Lateral Width	Weight
1	4170-0011-020	22 kg	27-31 cm	28-33 cm	40-48 cm	21-25 cm	19-28 cm	8.2 kg
2	4170-0021-060	32 kg	31-35 cm	34-39 cm	48-55 cm	25-29 cm	23-32 cm	8.9 kg
3	4170-0031-090	75 kg	36-40 cm	40-46 cm	55-64 cm	30-36 cm	29-40 cm	12.3 kg