



Lasse upright and supine stander – Sizes Mini to 2

REHATEC®

The height adjustable upright and supine stander Lasse moves young users from a supine to a standing position and enables them to stand at eye level with their friends. The variable transfer height ensures a safe, stress-free positioning of the user in supine.



Lasse size 1
Frame colour: pure white
Fabric: light blue/grey




Lasse size 2
height adjustable
Frame colour: anthracite grey
Fabric: apple green/grey



Degree indicator for height
adjustable frame



Leg positioning,
angle adjustable



Lateral supports,
flip away



Homepage



Product page/
technical data

 **Paragon**
mobility

GMFCS
Level
IV - V



3 Sizes



0.65 m -
1.35 m



Max.
60 kg



6 cover
colours



22 frame
colours

**Also available
with height adjustment**



Infants



Children



Teenagers





Lasse upright and supine stander – Sizes 3 to 5

REHATEC®

The recline of the Lasse upright and supine stander can be steplessly adjusted and gently brings users of all ages to a standing position. Since the upright and supine stander adapts to users from head to toe, Lasse perfectly meets individual needs.



Lasse size 3
Frame colour: matt black
Fabric: green/grey




Lasse size 4
Frame colour: matt black
Fabric: yellow/grey



Individual leg abduction



Height adjustable for
transfer & improved access



Knee cups, contoured



Homepage



Product page/
technical data

 **Paragon**
mobility

GMFCS
Level
IV - V



3 Sizes



1.28 m -
2.00 m



Max.
130 kg



10 cover
colours



22 frame
colours

**Also available
with height adjustment**



Teenagers



Adults

