

Little Wave "Arc"



AHEAD OF THE **Curve**

Little Wave **"Arc"**



GET OUT & **GO!**

Arc is the lightest pediatric folding tilt-in-space wheelchair available. It's intuitive, easy-to-fold design and compact wheelbase make its folded footprint very small. It's the perfect fit, leaving plenty of room for a family on the go.


TAKE IT **EASY!**

Keeping parents in mind, the Little Wave Arc is your best solution for kids who need tilt-in-space. Our **Arc Tilt System** (patent-pending) results in effortless repositioning and a shorter, more maneuverable, wheelbase.



READY, SET, **GROW!**

Unlike shoes, you can't just go out and get another pair when your child outgrows their chair. So, we built in plenty of growth and offer unique programs for growing kids. We also designed the Arc with the latest high-strength, lightweight materials to make it durable.



EXPRESS YOURSELF

Everyone has a style that's their own and we want to make sure you can express yours. That's why we offer so many style options, including 2 color frames for free. Pick your style and show it off!



OPTIONS & ACCESSORIES



Folding Rigidizer Bar



Single Hand Tilt




Attendant Wheellock



Folding Angle Adjustable Push Handle



Foot Tilt



Canopy

SPECIFICATIONS

Transport Weight	20.5 lbs*	9.3 kg*
Seat Width	10-18 in.	25-46 cm
Seat Depth	12-20 in.	30-51 cm
Seat Height	14.5- 19 in.	37-48 cm
Weight Capacity	165 lbs.	75 kg
WC-19 Transit	Yes	Yes

* 14" x 14" Frame with Attendant Wheel Locks. Armrests, Rear Wheels and Footrests Removed.



Stevens Point, WI
800-981-1540
www.kimobility.com
DCN0827.0