

CASTER GLIDE™ | INSPIRED DETAILS

Our focus on end-user performance is evident even in the smallest details. An improved 7" caster design and adjustment-free flutter damping device means you will move more confidently at low speeds, in tight areas, or even while using a variety of alternative drive controls. You can now turn the way you want with full caster clearance even at 90° legrest position. This level of precision engineering helps ensure improved drive performance and a more consistent driving experience for all users.

FREE TO MOVE  
7" PERFORMANCE CASTER  
DESIGN + ADJUSTMENT  
FREE FLUTTER DAMPING

COLOUR OPTIONS

 Radiant Red

 Power Pink


 Ocean Blue


 Sunburst Orange


 Galactic Green

 Midnight Black


HYDROGRAPHIC OPTIONS


 Mossy Oak® Breakup® Country™

 CDG Poly

 Carbon Fiber

FABRIC OPTIONS

 Leatherette

 Mesh

SPECIFICATIONS

Maximum User Weight	136 kg
Maximum Speed	10 km/h
Range	21 km*
Base Width	24"
Base Length	35.5"
Minimum Turning Radius	20.5"
Ground Clearance	3"
Recommended Battery Type	Group 34, Group 24
Weight Including Batteries	175 kg (Gr 34), 184 kg (Gr 24)
Drive Electronics	R-Net 120A
Seat-to-Floor Height	17.5", 18.5", 19.5"

Optional Seat Elevator	12"
Tilt Options	0° – 50° Posterior 5°, 10° & 20° Anterior
Recline Options	85° – 120° (manual) 85° – 180° (power)
Legrest Angle	85° – 170°
Backrest Height	20", 23" – 28" (by 1" incr.)
Seat Depth	14" – 22" (by 1" incr.)
Seat Width	17" – 23" (by 2" incr.)
Armrest Pad Lengths	10", 13", 16", 18"

*\*Specifications may vary based upon configuration.  
\*Actual driving range will vary based on driving and battery conditions.*

permobil

YOUR WORLD  
WITHIN REACH.

M3  
Corpus





ADDED STABILITY WITH UNIQUE  
SUSPENSION GEOMETRY AND  
LINEAR SPRING DESIGN

PATENT PENDING FLEXLINK™ SUSPENSION  
FOR ADVANCED VIBRATION REDUCTION

## FLEXLINK™ | FEEL THE DIFFERENCE

Conquer your day with a suspension meticulously engineered to help you maintain better positioning, feel more secure, and reduce driving fatigue. The unique suspension geometry and linear spring design of the M3 Corpus provides a more stable chassis while the vibration damping FlexLink™ ensures a smoother ride over a variety of surfaces.

## ACTIVEREACH™ | REACH FOR MORE

Your world isn't one dimensional. Your seat elevator shouldn't be either. ActiveReach™ technology is the ultimate advantage in accessing your environment. By combining 12" of seat elevation and up to 20° of forward tilt, Permobil's ActiveReach™ technology enables you to reach beyond obstructions, into high cabinets, and can even assist in transfers. Truly putting your world within reach.



## HIGHER STANDARDS

This innovative combination of swing arms, linear spring design and FlexLink™ suspension enables the M3 to achieve full suspension at ANY elevation, seating position or driving speed so that you never have to sacrifice comfort for function.

For even more stability in elevation, the seating system moves rearward 3.5" and maintains three points of contact with the chassis at all times.

5 KM/H @ 12" OF  
ELEVATION + FULL  
SUSPENSION AT ANY  
SPEED, POSITION, OR  
ELEVATION

12" OF ELEVATION + 20°  
FORWARD TILT = ULTIMATE  
ACCESSIBILITY.

