



# Noah hi-lo base

# REHATEC®

Infants



Children



Teenagers



The Noah hi-lo base is a tilting frame for all common seat shells. Noah offers gas spring, hydraulic or powered height adjustment for an optimal seat height at all times. With outdoor tyres, it is also suitable for outdoor activities.

GMFCS  
Level  
IV - V

2 Sizes



Max.  
75 kg



22 frame  
colours

Noah size 2  
Frame colour: matt black



Outdoor frame



Height adjustment to  
lower the seat down to the  
ground



Compatible with the  
TRIXI seating system



Homepage



Product page/  
technical data

