



# Raz-AT

Australian Order Forms  
Effective September 1, 2020

Date: \_\_\_\_\_

Dealer / Location: \_\_\_\_\_ T: \_\_\_\_\_ F: \_\_\_\_\_

PO# / Client: \_\_\_\_\_ Contact: \_\_\_\_\_ Email: \_\_\_\_\_

**Raz-AT Mobile Shower Commode Chair** with 40° tilt range, height-adjustable stainless steel frame, four dual-locking 5" casters, stainless steel head support fixture, molded small head support pad, flip-up padded arm supports, foot supports. Exclusive to Raz, all seats include the Ischial Pelvic Alignment System [IPAS], which allows for 2" fore/aft adjustment. **User weight capacity 400 lbs.** Please enter special set-up requests, such as arm support height or tilt range restriction, in Notes [Sec.15].

## 1. Chair Codes by size – Must Choose One

Z300S18	Raz-AT w/ 18"W seat and base frame – Standard size, fits over most toilets <sup>1</sup>	ADD
Z300S__*	Non-Standard Width Raz-AT – Must add upcharge for custom width: see below * add custom width to complete code (16" / 20" / 22" / 24") – <b>Adjusta-Back must be selected for all custom widths [Sec.5]</b>	ADD
Z3SW16	Upcharge for custom 16"W seat and base frame <sup>2</sup> – will not fit over most toilets	ADD
Z3SW20	Upcharge for custom 20"W seat and base frame <sup>3</sup>	ADD
Z3SW22	Upcharge for custom 22"W seat and base frame <sup>4</sup>	ADD
Z3SW24	Upcharge for custom 24"W seat and base frame <sup>4</sup>	ADD

<sup>1</sup> Fits all 18"W seats

<sup>2</sup> Visco-Foam Interface seats can be modified to fit.

<sup>3</sup> 18"W seats can be selected but must also select Seat Support Adapter Kit Z4SSKT [Sec.2]

<sup>4</sup> Z4V23N Visco Foam Interface HD seat can be modified to fit. Seat Support Kit (Z4SSKT) required for 24"W seat frame.

## Custom Seat / Base Frame Length and Seat Height

Z3SFL3	Upcharge for custom 3" Longer Seat Frame: required for seats ≥ 21"D	ADD
Z3BFL3	Upcharge for custom 3" Longer Base Frame: required for seats ≥ 22"D (typically used in conjunction with Z3SFL3)	ADD
Z3SL5	Upcharge for custom Super-Low Frame (5" lower) – height adjustment range (molded seat): 16"-20" <sup>1</sup>	ADD

<sup>1</sup> To lower height adjustment range by 1" (15"-19"), select 4" caster upgrade option Z64T (will not fit over most toilets)

## 2. Seats – Standard

Molded Seats (thickness ~2.5"-3")			Visco-Foam Interface Seats (thickness ~3.5")		
Z405	18"W x 16"D – Front Access	ADD	Z436BRI	18"W x 16"D – Front Bridged	ADD
Z410	18"W x 16"D – Front Bridged	ADD	Z436FRO	18"W x 16"D – Front Access	ADD
Z415	18"W x 19"D – Front Access	ADD	Z436LTO	18"W x 16"D – Left Access	ADD
Z420	18"W x 19"D – Front Bridged	ADD	Z436RTO	18"W x 16"D – Right Access	ADD
<b>Specialty Seats</b>			Z436RRO	18"W x 16"D – Rear Access	ADD
Z4SOM	Shower-Only Seat – 18"W x 16"-21"D	ADD	Z439BRI	18"W x 19"D – Front Bridged	ADD
Z4EJRM	E&J Replica Seat – Molded – 18"W x 17"D	ADD	Z439FRO	18"W x 19"D – Front Access	ADD
Z4EJRV	E&J Replica Seat – Visco Foam – 18"W x 17"D	ADD	Z439LTO	18"W x 19"D – Left Access	ADD
Z4V23N	23"W x 18"D HD – Front Access Notched at rear to maintain fore/aft adjustment range (IPAS) – 22"W frame	ADD	Z439RTO	18"W x 19"D – Right Access	ADD
			Z439RRO	18"W x 19"D – Rear Access	ADD
			Z4SSKT	Seat Support Kit for Raz-AT/CAT (to support seats that are narrower than seat frames)	ADD

### 3. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required.

Z4CS	Custom Seat <b>Select if both width and depth are <math>\leq 19''</math></b>	ADD
Z4CSHD	Custom Seat <b>Select if either width or depth are between <math>20'' \leq 24''</math></b>	ADD
Z4CSHDX	Heavy Duty Custom Seats <b>Select if width and/or depth is <math>\geq 24''</math></b>	ADD

#### MUST FILL OUT PARAGON CUSTOM SEAT FORM

- Custom Seat width should not be narrower than frame width.
- Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support without ordering wider chair frame.
- The Ischial Pelvic Alignment System (IPAS) allows 2" fore/aft adjustment of the seat (seat depth and aperture position) on all seats.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, etc.

### 4. Commode Set up – These set up configurations are for dealer use only. Commodes will not be supplied out of the box with the below measurements.

**Set armrest height to:** \_\_\_\_ Measured from underside of the seat to the top of armrest.

**Seat toilet clearance to:** \_\_\_\_ This is the minimum clearance between the rear bar of the commode and the ground. (Excludes Zimpan Kit)

**Footrest length:** \_\_\_\_ Measured from underside of the seat to the top of footplate.

### 5. Adjusta-Back Frames. Where your back frame is different from seat frame, please also refer to [Sec.6]




Includes height- and angle-adjustable arms, tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (85°/90°/95°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. 24"H Adjusta-Back includes straight posts, on which the upholstery cannot be raised.

Standard Width Adjusta-Backs (18"W)			Non-Standard Width Adjusta-Backs		
Z50918	Adjusta-Back – 15"H	ADD	<b>* Complete Adjusta-Back order code by adding desired back frame width (16"/20"/22"/24") Note: 18" is the standard width.</b>		
Z51018	Adjusta-Back – 17"H	STD	Z509__*	Adjusta-Back – 15"H non-std width	ADD
Z51518	Adjusta-Back – 20"H	ADD	Z510__*	Adjusta-Back – 17"H non-std width	ADD
Z52018	Adjusta-Back – 24"H	ADD	Z515__*	Adjusta-Back – 20"H non-std width	ADD
Z5DTL	Dual Triggers – both on left side <sup>1</sup>	N/C	Z520__*	Adjusta-Back – 24"H non-std width	ADD
Z5DTR	Dual Triggers – both on right side <sup>1</sup>	N/C	Z515ATR__*	Reclining Adjusta-Back <sup>2</sup> 30° recline, 20" back height with 2-Point Reclining Arms (tilt is limited to 25°)	ADD
Z515ATR18	Reclining Adjusta-Back <sup>2</sup> 30° recline, 20" back height with 2-Point Reclining Arms (tilt is limited to 25°)	ADD	Z520ATR__*	Reclining Adjusta-Back <sup>2</sup> 30° recline, 24" back height with 2-Point Reclining Arms (tilt is limited to 25°)	ADD
Z520ATR18	Reclining Adjusta-Back <sup>2</sup> 30° recline, 24" back height with 2-Point Reclining Arms (tilt is limited to 25°)	ADD			

<sup>1</sup> Dual Tilt Triggers option positions both triggers on one side, allowing for single-handed activation. Available for all back widths.

<sup>2</sup> On reclining backs (Z515ATR/Z520ATR), Dual Triggers for recline are mounted on the left side and for tilt on the right side.

**6. Custom Back Frame Widths Different From Seat Frames – Must choose back frame [Sec.5] (configurations below are not available for recliner backs)**

Z5NB__ __ *	Custom Back – 2” Narrower than Seat Frame	ADD	Seat Frame Width Back Frame Width	16” <sup>1</sup> 14”	18” <sup>2</sup> 16”	20” <sup>3</sup> 18”	22” <sup>4</sup> 20”	24” <sup>5</sup> 22”	
* Enter seat and back frame widths to complete code. eg. Z5NB2018 specifies a 20”W seat frame with an 18”W back frame									
Z5NB__ __ *	Custom Back – 4” Narrower than Seat Frame – requires custom seat	ADD	Seat Frame Width Back Frame Width	N/A N/A	18” <sup>1</sup> 14”	20” <sup>3</sup> 16”	22” <sup>4</sup> 18”	24” <sup>5</sup> 20”	
* Enter seat and back frame widths to complete code. eg. Z5NB2016 specifies a 20”W seat frame with a 16”W back frame									
Z5WB__ __ *	Custom Back – 2” Wider than Seat Frame	ADD	Seat Frame Width Back Frame Width	16” 18”	18” 20”	20” <sup>3</sup> 22”	22” <sup>6</sup> 24”	N/A N/A	
* Enter seat and back frame widths to complete code. eg. Z5WB1820 specifies an 18”W seat frame with a 20”W back frame									
Consider also using Arm Spacer Kit (Z718/Z719) as another solution for providing more width between the arms.									
<sup>1</sup> Requires Custom Seat <sup>2</sup> Compatible with all 18”W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats. <sup>3</sup> Custom Seat recommended. Also compatible with all 18”W seats, requires Seat Support Kit Z4SSKT [Sec.2] <sup>4</sup> Custom Seat recommended. Also compatible with modified Z4V23N (18” seat depth with limited adjustability). <sup>5</sup> Custom Seat recommended. Compatible with modified Z4V23N (18” seat depth with limited adjustability, requires Seat Support Kit (Z4SSKT) [Sec.2] and modified Z4V26N (18” seat depth with limited adjustability). <sup>6</sup> Custom Seat recommended. Also compatible with Z4V23N (18” seat depth).									

**7. Back Support Options – These back supports replace the standard, tension-adjustable upholstery and are recommended for infection-control environments.**

Z5S14	14”W Symphony Back <sup>1</sup> – recessed solid back	ADD	Z749	Harmony Back <sup>2</sup> (contoured solid back – 18”W x 13”H) Fits 17” or 20” back heights, 18”W back frames only	ADD
Z5S16	16”W Symphony Back <sup>1</sup> – recessed solid back	ADD			
Z5S18	18”W Symphony Back <sup>1</sup> – recessed solid back	ADD			
Z5S20	20”W Symphony Back <sup>1</sup> – recessed solid back	ADD			
Z5S22	22”W Symphony Back <sup>1</sup> – recessed solid back	ADD			
Z5S24	24”W Symphony Back <sup>1</sup> – recessed solid back	ADD			
Z5S26	26”W Symphony Back <sup>1</sup> – recessed solid back	ADD			
<sup>1</sup> Not compatible with 15”H back – all Symphony Back panels are 13”H <sup>2</sup> Not recommended for reclining backs.					

## 8. Arm Supports and Accessories

Z5BAR13	Flip-Up Arm Supports – 13”L pad (pr)	STD	Z7PAML	Pivoting Arm Mount – Left (for Z7FL/Z7TR)	ADD
Z710	Armrest Lock Kit (for Z505/Z506) (pr)	ADD	Z7PAMR	Pivoting Arm Mount – Right (for Z7FL/Z7TR)	ADD
Z52PAT	2-Point Flip-up Arm Supports <sup>1</sup> (pr) (standard with reclining backs)	ADD	Z718	Arm Support Spacer Kit – 1/2” per side (pr)(for Z509/Z510/Z515/Z520)	ADD
Z7FLL	Long Molded Flat Armpad <sup>2</sup> – Left (14”L x 3.25”W)	ADD	Z719	Arm Support Spacer Kit – 1” per side (pr)(for Z509/Z510/Z515/Z520)	ADD
Z7FLR	Long Molded Flat Armpad <sup>2</sup> – Right (14”L x 3.25”W)	ADD	Z7PHG	Pivoting Hand Grips (pr)	ADD
Z7TRL	Molded Arm Trough <sup>2</sup> – Left (18.5”L x 4.5”W)	ADD	Z7APS_ _*	Anterior Posture Support Bar * Enter back frame width to complete code (14”/16”/18”) <sup>3,4</sup>	ADD
Z7TRR	Molded Arm Trough <sup>2</sup> – Right (18.5”L x 4.5”W)	ADD	Z7APS_ _*	Anterior Posture Support Bar * Enter back frame width to complete code (20”/22”/24”) <sup>3,4</sup>	ADD
<sup>1</sup> For chairs with same back width and frame width only <sup>2</sup> May require Z7PAM for arm support to flip up and/or for lateral supports to swing away			<sup>3</sup> Not available with Z7FL / Z7TR <sup>4</sup> Special length required with Z718 (+1”) and Z719 (+2”).		

**9. Foot and Leg Supports** – Standard Foot Supports (Z705N) are included with the chair. The selection of any Foot / Leg Support Option below automatically replaces the relevant Standard Foot Support. e.g. By selecting Elevating Leg Support – Left (Z7ELRL), the chair will be supplied with one Foot Support (Z705N) on the right side and a Z7ELRL on the left side.

Z705N	Foot Supports 16”-20.5” (pr)	STD	Z7VS	V-Style Foot Support <sup>3</sup> (1-piece)	ADD
Z705N16	Foot Supports for 16”W frames 16”-20.5” (pr)	N/C	Z7ELRL	Elevating Leg Support <sup>4</sup> – Left (includes Z7ADM)	ADD
Z7FET2	Footrest Extension Tubes – 2” longer <sup>1</sup> (pr)	ADD	Z7ELRR	Elevating Leg Support <sup>4</sup> – Right (includes Z7ADM)	ADD
Z7MFX8	MFX8 Foot Supports 8”-12” (pr)	ADD	Z7ELROFL	Forward Offset Elevating Leg Support <sup>5</sup> – Left (includes Z7ADM)	ADD
Z7MFX12	MFX12 Foot Supports 12”-16.5” (pr)	ADD	Z7ELROFR	Forward Offset Elevating Leg Support <sup>5</sup> – Right (includes Z7ADM)	ADD
Z7ADM	Angle/Depth Adj Molded Footplates <sup>2</sup> (pr)	ADD	Z7RLSL	Universal Residual Limb Support – Left	ADD
Z7FET2A	Footrest Extension Tubes for Z7ADM (pr) (2” longer)	ADD	Z7RLSR	Universal Residual Limb Support – Right	ADD
Z7ADABS	Angle/Depth Adj ABS Footplates (pr) For 16”W frames and Z7MFX8	ADD	<sup>3</sup> Custom required for sizes other than 18”W seat frames <sup>4</sup> Z7ELRL/R compatible with frame widths ≥ 18” and adjusted seat depths ≤ 18” or ≤ 21” w/ 3” longer seat frames (Z3SFL3) <sup>5</sup> Z7ELROFL/R compatible with frame widths ≥ 18” for adjusted seat depths 19”-20” on standard frame or 22”-23” w/ 3” longer seat frames (Z3SFL3)		
Z7FBT	Flip-Back Foot Support (1-piece) with calf pads <sup>3</sup>	ADD			
<sup>1</sup> Not compatible with Z7MFX8 / Z7ADM / Z7VS / Z7FBT <sup>2</sup> Not compatible with Z300S16 or Z7MFX8					

**10. Foot and Leg Support Accessories**

Z7FSLK	Foot Support Lock (pr)	ADD	Z708ADJ	Adjustable Heel Loops (pr) – not available for Z7VS / Z7FBT	ADD
Z7CS_*	Snap-On Calf Strap – Bodypoint Aeromesh	ADD	Z7NCS	Neoprene Covers for Standard Footplates (pr)	ADD
	* Enter seat frame width to complete code (16"/18"/20"/22"/24")		Z7NC16	Neoprene Covers for narrow molded footplates (pr)	ADD
Z737CPL	Calf Panel – Bodypoint Aeromesh <sup>1</sup> Fits frame widths 16"/18"/20"	ADD	Z7NCA	Neoprene Covers for Angle/Depth Adj Molded Footplates (pr)	ADD
Z706	H-Strap – for 18"/20"W frames – hook and loop attachment	ADD	Z7NCABS	Neoprene Covers for ABS Adj Angle/Depth Footplates Z7ADABS (pr)	ADD
Z712	H-Strap – for 22"/24"W frames – hook and loop attachment	ADD	Z7OFLM	Lateral Offset Foot Support Receiver – moves foot support 2" laterally or medially Quantity (1 or 2) ___	\$160
<sup>1</sup> Not compatible with Z7MFX8 / Z7FBT / Z7VS					

**11. Positioning Pads**

Z7HSTS	Head Support Pad – Small (10" x 4")	STD	Z7LGT1L	Leg Adductor w/ Small Pad <sup>1</sup> – Left	ADD
Z7HSTL	Head Support Pad – Large (12" x 6")	ADD	Z7LGT1R	Leg Adductor w/ Small Pad <sup>1</sup> – Right	ADD
Z7OHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment	ADD	Z7LGT2L	Leg Adductor w/ Large Pad <sup>1</sup> – Left	ADD
Z7LTS1L	Lateral Thoracic Support w/ Small Pad <sup>1,2</sup> – Left	ADD	Z7LGT2R	Leg Adductor w/ Large Pad <sup>1</sup> – Right	ADD
Z7LTS1R	Lateral Thoracic Support w/ Small Pad <sup>1,2</sup> – Right	ADD	Z7HP1L	Hip Adductor w/ Small Pad <sup>1,3</sup> – Left	ADD
Z7LTS2L	Lateral Thoracic Support w/ Large Pad <sup>1,2</sup> – Left	ADD	Z7HP1R	Hip Adductor w/ Small Pad <sup>1,3</sup> – Right	ADD
Z7LTS2R	Lateral Thoracic Support w/ Large Pad <sup>1,2</sup> – Right	ADD	Z7HP2L	Hip Adductor w/ Large Pad <sup>1,3</sup> – Left	ADD
Z742L	Lateral/Adductor Extension Plate (5.5" x 3.5") – Left QTY ___	ADD	Z7HP2R	Hip Adductor w/ Large Pad <sup>1,3</sup> – Right	ADD
Z742R	Lateral/Adductor Extension Plate (5.5" x 3.5") – Right QTY ___	ADD	Z7MTST1	Medial Thigh Support – Small	ADD
			Z7MTST2	Medial Thigh Support – Large	ADD
<sup>1</sup> Adjustable Pad Sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3".					
<sup>2</sup> Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For additional medial offset, specify Z742 extension plate.					
<sup>3</sup> Fits Adjusta-Back only					

## 12. Pelvic and Chest Belts – A chest belt should not be used without a pelvic belt

Pelvic Belts				Chest Belts					
ZPBS	Pelvic Belt – Small (frame widths 14"-16", 2-piece w/ side-release buckle)	ADD		ZCBS	Chest Belt <sup>2</sup> – Small (back widths 14"-16", 1-piece w/ hook and loop closure)	ADD			
ZPBM	Pelvic Belt – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	ADD		ZCBM	Chest Belt <sup>2</sup> – Medium (back widths 16"-18", 1-piece w/ hook and loop closure)	ADD			
ZPBL	Pelvic Belt – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	ADD		ZCBL	Chest Belt <sup>2</sup> – Large (back widths 18"-22", 1-piece w/ hook and loop closure)	ADD			
ZPBXL	Pelvic Belt – XLarge (seat frame widths 22"-30", 2-piece w/ side-release buckle)	ADD		ZCBXL	Chest Belt – XLarge <sup>2</sup> (back frame widths 22"-30", 1-piece w/ hook and loop closure)	ADD			
ZPBICS	Infection Control Pelvic Belt <sup>1</sup> – Small (frame widths 14"-16", 2-piece w/ side-release buckle)	ADD		ZCBICS	Infection Control Chest Belt <sup>1,2</sup> – Small (back widths 14-16", 1-piece w/ side-release buckle)	ADD			
ZPBICM	Infection Control Pelvic Belt <sup>1</sup> – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	ADD		ZCBICM	Infection Control Chest Belt <sup>1,2</sup> – Medium (back widths 16"-18", 1-piece w/ side-release buckle)	ADD			
ZPBICL	Infection Control Pelvic Belt <sup>1</sup> – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	ADD		ZCBICL	Infection Control Chest Belt <sup>1,2</sup> – Large (back widths 18"-22", 1-piece w/ side-release buckle)	ADD			
ZPBICXL	Infection Control Pelvic Belt <sup>1</sup> – XLarge (seat frame widths 22"-30", 2-piece w/ side-release buckle)	ADD		ZCBICXL	Infection Control Chest Belt <sup>1,2</sup> – XLarge (back frame widths 22"-30", 1-piece w/ side-release buckle)	ADD			
ZPBBPS	Pelvic Belt – Bodypoint Aeromesh – Small	ADD		ZCBBPS	Chest Belt – Bodypoint Aeromesh – Small	ADD			
ZPBBP_*	Pelvic Belt – Bodypoint Aeromesh * enter size to complete code (M / L)	ADD		ZCBBP_*	Chest Belt – Bodypoint Aeromesh * enter size to complete code (M / L / XL)	ADD			
ZPBBPXL	Pelvic Belt – Bodypoint Aeromesh – XLarge	ADD							
<b>Recommended sizes for Bodypoint belts</b>	Chest belt	ZCBBPS	14"-16"W back frame	ZCBBPM	16"-18"W back frame	ZCBBPL	18"-22"W back frame	ZCBBPXL	22"-30"W back frame
	Pelvic belt	ZPBBPS	14"-16"W seat frame	ZPBBPM	16"-18"W seat frame	ZPBBPL	18"-22"W seat frame	ZPBBPXL	22"-30"W seat frame
<sup>1</sup> Non-absorbant belts – recommended for infection-control environments									
<sup>2</sup> Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.									

**13. Other**

Z675	Directional Caster Lock – fits on one rear caster (facilitates steering)	ADD	ZIHW	Whizard – winged urine deflector <sup>2</sup>	ADD
Z66T4	Anti-Tippers for 4" casters – Rear (pr)	ADD	ZBUDBH__*	Bard Urine Drainage Bag Holder * Enter size to complete code (16"/18"/20"/22"/24")	ADD
Z66T5	Anti-Tippers for 5" casters – Rear (pr)	ADD			
Z66T6	Anti-Tippers for 6" casters – Rear (pr)	ADD	ZSPKT22	AT Self-Propel Kit with 22" wheels – incl anti-tippers (tilt limited to 25°)	ADD
Z64T	4" Caster Upgrade (4x)	ADD			
Z66T	6" Caster Upgrade (4x)	ADD	ZSPKT24	AT Self-Propel Kit with 24" wheels – incl anti-tippers (tilt limited to 25°)	ADD
ZPSK1	Commode Pan Spacer Kit <sup>1</sup> – 1" lower	ADD	ZIMPANKIT	Commode pan and seat attachments	ADD
ZPSK2	Commode Pan Spacer Kit <sup>1</sup> – 2" lower	ADD			
ZPSK3	Commode Pan Spacer Kit <sup>1</sup> – 3" lower	ADD	ZSCLAMPKIT	Seat Clamp Kit. (Only required if not ordering a seat from Paragon)	ADD
ZPSK4	Commode Pan Spacer Kit <sup>1</sup> – 4" lower	ADD			
ZTLS	Tilt Limit Spacers – 1/2" spacers snap onto gas struts to limit tilt (each reduces tilt range approx. 5°. Can be stacked for additional tilt limiting. Quantity ____	N/C			

<sup>1</sup> Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Reduces toilet bowl clearance height.

<sup>2</sup> Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.

**Some assembly of Raz products is required.**

**14. Notes:**

## Specifications

## Raz-AT (Attendant Tilt)

Weight Capacity	400 lbs
Overall Width	Seat width plus 4.5" with Standard Arm Pads
Overall Length (with foot supports and headrest adjusted to middle positions)	
- at 0° tilt and standard foot supports	43"
- at 40°	54"
Seat Depth	16" or 21" depending on seat/cushion length
Toilet Clearance Height (under seat)	
- with pan brackets	17"-21" (1" increments)
- without pan brackets	17.75"-21.75 (1" increments)
Toilet Clearance Width	15.75" for standard 18"W seat frame. For other widths, seat frame width minus 2.25"
Arm Support Height - above molded seat (1/2" increments)	
- 17"H Fixed Back	8.5"
- <b>20"H Fixed Back - standard</b>	<b>9.5"</b>
- 15"H Adjusta-Back	5.5"-8.5"
- 17"H Adjusta-Back	6.5"-9.5" (1/2" increments)
- 20"H Adjusta-Back	7.5"-10.5" (1/2" increments)
- 24"H Adjusta-Back	7.5"-10.5" (1/2" increments)
- 20"H Reclining Adjusta-Back	7.5"-10.5" (1/2" increments)
- 24"H Reclining Adjusta-Back	7.5"-10.5" (1/2" increments)
Seat-to-Floor Height (with molded seat)	21"-25" in 1" increments
Product Weight (with 18"W seat)	61 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)