



Raz-AT600

Attendant Propel Heavy Duty Tilt
Rehab Shower Commode 600-lb Capacity

Raz-AT600 [Attendant Tilt, 600-lb (275 kg)] Heavy Duty Rehab Shower Commode Chair with 25° tilt range. The specially engineered, Weight Assist Spring System (WASS) is available in two weight ranges and is adjustable to accommodate the user's weight and center of mass for tilt operation with minimal effort by the attendant.



FEATURES

Extensive range of options available

23"W (58 cm) x 18"D (46 cm) contoured, molded seat

Exclusive to Raz, Ischial Pelvic Alignment System (IPAS), provides 2" (5 cm) of fore/aft or 1" (2.5 cm) of lateral seat adjustment

22"W (56 cm) stainless steel frame with 4" (10 cm) of height adjustment

Adjustable heavy-duty footrests

Dual-locking (swivel and roll) 6" (150 mm) casters



Tension-adjustable back upholstery

Flip-up arms

Adjustable back angle

Backrest assembly adjusts fore/aft for seat depth change and center-of-mass adjustment

25° of tilt

ALSO AVAILABLE

- Frame widths up to 36" (91 cm)
- Custom frame depths
- Custom seats

SPECIFICATIONS	Raz-AT600 (Attendant Tilt 600 lbs capacity)
WEIGHT CAPACITY	600 lbs (275 kg)
FRAME WIDTH (seat and back frame)	22" (56 cm)
SEAT WIDTH (commode seat)	22" - 30" (56 cm - 76 cm)
SEAT DEPTH (front of back post to front of seat)	18" (46 cm)
OVERALL WIDTH (with standard armpad)	Seat frame width + 4.5" (11 cm)
OVERALL LENGTH With Z709 Footrests	38" - 39" (96 - 99 cm)
CLEARANCE HEIGHT UNDER SEAT	with pan brackets: 18" - 19" (43 - 48.5 cm) without pan brackets: 18.75" - 19.75" (47.5 - 50 cm)
ARMREST HEIGHT ABOVE SEAT	9.5" (24 cm)
SEAT-TO-FLOOR HEIGHT	21.75" and 22.75" (55 cm and 58 cm)
PRODUCT WEIGHT	5" casters: 82 lbs (37 kg) 6" casters: 84 lbs (38 kg)
FRAME AND HARDWARE MATERIAL	Stainless Steel (Lifetime warranty against corrosion)

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions.

Specifications, and options are subject to change without prior notice.

MANUFACTURED BY:

RAZ Design Inc. 22 Howden Road, Toronto, Ontario Canada,

M1R 3C7 T: (877) 720-5678 E: info@razdesigninc.com

RAZSALGO05r01 Rev. Date: 2017-02-27

Raz Design products are distributed in Australia by:

Paragon Mobility

3/10 Tapnor Crescent, Brendale, QLD, 4500 Australia

Tel: 13000 652 382

Fax: (07) 3112 6125

www.paragonmobility.com.au

