



Zwim

Pool Chair

Zwim offers a safe way to enter and exit a pool that has a ramp or for transportation from the locker room to a pool lift. A contoured seat (with drainage slot) provides comfort. Ideal for chlorinated environments.



330
lbs

Weight Capacity

STANDARD FEATURES



SPECS

SPECIFICATIONS	Zwim POOL CHAIR
Weight Capacity	330 lbs (150 kg)
Seat Width	19.5" (50 cm)
Seat Depth (front of back post to front of seat)	17" (43 cm)
Back Width	18.5" (47 cm)
Back Height (from top of seat to top of upholstery)	17" (43 cm)
Arm Support Height Above Seat	9" (23 cm)
Seat-To-Floor Height* (2.5 cm [1"] increments)	19" - 23" (48 - 58 cm)
Overall Width	26.5" (67 cm)
Product Weight	39.5 lbs (18 kg)
Frame Material	Aluminum

Raz Design products are distributed in Australia by:
Paragon Mobility

3/10 Tapnor Crescent, Brendale, QLD, 4500 Australia
 Tel: 13000 652 382
 Fax: (07) 3112 6125
www.paragonmobility.com.au

