

ROGUE



YOU'RE YOUNG.  
YOU'RE ACTIVE.  
YOU'RE YOU.



GROWTH WITH AN EDGE

## ROGUE XP

Rogue XP is the ultimate performance chair; encouraging you to embrace your personal style at the age when you care the most. It gives you the necessary built-in growth; no need to compromise fit, function, ride or style.

*Focus your effort on pushing your limits... not your chair.*

Whether you're climbing into the driver's seat or hitching a ride with friends, Rogue XP's **LOCKING BACK MECHANISM** and **11.6 LB. TRANSPORT WEIGHT** makes it easier.



## CREATIVE FREEDOM IS YOURS

Rogue XP is your blank canvas giving you the freedom to choose your frame, accent, release cable and wheel color. Whether you're unique, chic, edgy or minimalistic... create your ID with Rogue XP.

## GROWTH IS YOURS

Optimize fit, function and ride with up to 3" of depth and width growth. Rogue XP changes with you over time so outgrowing your chair is a problem of the past.

## IDEAL POSITION IS YOURS

The integrated swing-away arm and backrest system adjusts independently allowing you to achieve ideal backrest and arm positioning.



## SPECS

Transport Weight*	11.6 lbs.	5.26 kg
Seat Width	10 - 18"	25.4 - 45.72 cm
Seat Depth	10 - 20"	25.4 - 50.8 cm
Front Seat Height	13.5 - 20.5"	34.29 - 52.07 cm
Rear Seat Height	11.5 - 20"	29.21 - 50.8 cm
Camber	0°, 2°, 6°, 8°	0°, 2°, 6°, 8°
Tubing Size	1 3/8"	3.49 cm
Weight Capacity	200 lbs.	90.72 kg
Ansi Resna WC-19 Transit Approved	Yes	Yes

\* 14x14" Standard Length Frame, Adjustable Height Back - 12", 4" Poly Casters, Tubular Open Footplate and Standard Upholstery. Rear Wheels Removed.