

Starlight Car Seat

Fully adjustable special purpose car restraint, offering excellent posture support options and an inbuilt tilt function. Starlight is suitable for young children from 9kg, and will grow with the child right up to 36kg.

Excellent postural support

Inbuilt adjustment to the rear of the seat allows the backrest height to be increased for growth

Large abdominal pad **protects child** at harness buckle

Crash tested to ECE R44-04



Uses **vehicle seatbelt** to secure the child

5 point comfort harness with shoulder and abdominal padding, and centre pull adjuster

ACCESSORIES

Forearm Support Tray
4141-5530-000



Abduction Pommel
4141-1510-000



Inbuilt Cup Holder

Vehicle Seating

Package includes

- Adjustable headrest
- Backrest infill padding
- Slim lateral supports
- 20° tilt
- Base board accessory
- Swivel base system - (Left turning) to enable turning of the car seat towards the doorway to help with transfers
- ISOFIX seat attachment
- Grey/Black upholstery



SLIM SWIVEL BASE

Allows more space and headroom inside the car when turning out.



Left Turning Swivel Package: 4141-0011-090

Code	Seat Width (pelvis)	Seat Width (knee)	Seat Depth (min-max)	Back Height	Crown Height	Height Laterals	Distance Laterals	Car Seat Weight	Capacity
4141-0011-090	30cm	35cm	26cm	27-39cm	51-62cm	17-27cm	14-19cm	8.5kg	9-36kg

*Reducer cushion available