

MADE IN GERMANY

max speed 10 km/h user capacity 140kg

range up to 20km

battery 6.2Ah, 36v, 224Wh motor 250W / 450W peak

12Nm / 25Nm peak

unit weight 7.2kg

- swiveling drive wheel enables all terrain use
- rigid and folding frame compatibility
- stop with mounted button
- stop with pushrims at low indoor settings
- adjust speed without disengaging
- easy attachment and removal
- USB C output to charge devices from smoov
- connect to smart phone with free app
- safe for air travel according to IATA



www.smoovbyalber.com.au