

TRACKZ

M O B I L I T Y

OFF THE BEATEN TRACKZ



We manufacture and sell the best wheelchairs for all your outdoor activities.

WE ARE THE
OFFROAD
SPECIALISTS!

TRACKZ

M O B I L I T Y

info@trackzmobility.com

  [TRACKZMOBILITY.COM](https://www.trackzmobility.com)

SPECS SHEET HP1

FRAME :	ALUMINUM
SEAT WIDTH :	14" TO 18"
SEAT DEPTH :	14" TO 18"
FRONT SEAT HEIGHT :	16" TO 20"
REAR SEAT HEIGHT :	14.5" TO 18.5"
SEAT TO FOOT REST :	12" TO 19"
SEAT BACK HEIGHT :	10.5" TO 15"
SEAT BACK ANGLE :	80 TO 105 DEG
CENTER OF GRAVITY :	1.5" TO 5"
COLOR ANODIZING PACKAGE :	NATURAL, BLACK, RED, PURPLE
FRONT WHEEL :	8" X 2"
REAR WHEELS :	24"
WHEELS LOCKS :	YES
BACK UPHOLSTERY COLOR :	BLACK, GREEN CAMO, PINK CAMO
ANTI-TIP :	OPT
PUSH HANDLES :	OPT
POWER ASSIST :	OPT



*WEIGHT: 21LBS *WITHOUT WHEELS AND ACCESSORIES



GET OUTSIDE AND **PLAY!**

TRACKZ
M O B I L I T Y

info@trackzmobility.com

  [TRACKZMOBILITY.COM](https://www.trackzmobility.com)