



Trakz

Tub Sliders

Trakz is a tub slider that consists of a mobility base, a tub base and a bridge. This system allows the user to be showered over a bathtub.



Trakz-AT (tilting)

FEATURES

Patented Zbrilok™ Safety System – Trakz is the world's first tub slider system with a patented safety feature that prevents the seat and user from moving onto the bridge unless the bridge is securely latched to both the tub base and the mobility base. This patented design, available only from Raz, can prevent serious injuries.



Rotational Tilt – Dual-Arc Tilt Geometry (Trakz-AT only)

- 35-degree tilt range – for pressure redistribution, trunk positioning and raising client's feet high enough to clear the tub wall
- Rotational tilt geometry reduces tilt effort and allows the Trakz to be mounted much closer to the rear wall

SPECIFICATIONS	Trakz-AP	Trakz-AT
WEIGHT CAPACITY		352 lbs (160 kg)
OVERALL WIDTH – with arm supports		23" (58.5cm)
– with arm supports removed		20.75" (53cm)
OVERALL LENGTH with standard footrests	35.5" (90.25cm)	36" (91.5cm)
TILT RANGE	N / A	35°
BACK SUPPORT HEIGHT		20" (51cm)
BRIDGE LENGTH	7", 8.5", 9.5", 10.5", 12.5", 14.5", 16.5" + custom lengths	
	18cm, 21.5cm, 24cm, 26.5cm, 32cm, 37.5cm, 42cm + custom lengths	
UNDERSIDE OF BRIDGE TO TOP OF SEAT	6.25" (15.75cm)	8.5" (21.5cm)
UNDERSIDE OF BRIDGE TO FLOOR		
– base configuration – high range	15.75"- 19.75" (40-50cm) in 1" (2.5cm) incr.	
– base configuration – low range	11.75"-15.75" (30-40cm) in 1" (2.5cm) incr.	
SEAT WIDTH		18" (46cm)
SEAT DEPTH		17" (43cm)
CLEARANCE HEIGHT UNDER SEAT		
– base configuration – high range	16.25"-20.25" (41-51.5cm) in 1" (2.5cm) incr.*	
– base configuration – low range	11.25"-15".25 (28.5-38.75cm) in 1" (2.5cm) incr.*	
* without brackets +0.75"		
ARMREST HEIGHT ABOVE SEAT		10" (25cm)
SEAT-TO-FLOOR HEIGHT		
– base configuration – high range	22"-25" (56-63.5cm) in 1" (2.5cm) incr.	24.5"-27.5" (62-70cm) in 1" (2.5cm) incr.
– base configuration – low range	17"-21" (43.25-53.5cm) in 1" (2.5cm) incr.	19.5"-23.5" (49.5-60cm) in 1" (2.5cm) incr.
PRODUCT WEIGHT	77.5 lbs (35kg)	95 lbs (43kg)
FRAME MATERIALS		
– base frame		Aluminum ¹
– back frame and footrests		Stainless Steel ¹
– hardware		Stainless Steel ¹

¹ Lifetime warranty against corrosion

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. Specifications and options are subject to change without notice.