

permobil
POWER WHEELCHAIRS

TOUGH ENOUGH



M300
CORPUS **HD**



RUGGED STYLE, REFINED DETAILS

M300 HD takes all the things you love about the Permobil's Mid Wheel Drive and amplifies them. Outfitted with stronger actuators, more powerful motors, and a reinforced seating system; then completed with elegant features like fully adjustable leg rests and up to four power seat functions, the M300 HD is tough enough.

COLOR OPTIONS



Volcano
Red



Popstar
Pink



Cobalt
Blue



Platinum
Silver



Sublime
Green



Purple
Rain



Phantom
Grey



Onyx
Black

FABRIC OPTIONS



Leatherette



Mesh

SPECIFICATIONS

| | |
|----------------------------|--|
| Maximum User Weight | 450 lbs |
| Maximum Speed | 5 mph |
| Range | 17.8 mi* |
| Base Width | 26" |
| Base Length | 38" |
| Minimum Turning Radius | 22" |
| Ground Clearance | 3" |
| Recommended Battery Type | Group 24 & Group 34 |
| Weight Including Batteries | 445 lbs with Group 24 425 lbs with Group 34 |
| Drive Electronics | R-Net 120A |
| Maximum Slope | 7.5° Up, 7.5° Down† |

| | |
|------------------------|---|
| Seat-to-Floor Height | 18", 19", 20" |
| Optional Seat Elevator | 8" |
| Power Tilt | 0° - 45° (power) |
| Recline Options | 85° - 120° (manual) 85° - 150° (power) |
| Legrest Elevation | 85° - 175° |
| Backrest Height | 21", 24" - 29" (by 1" incr.) |
| Seat Depth | 19" - 25" (by 1" incr.) |
| Seat Width | 19", 21", 23", 25" |
| Armrest Pad Lengths | 13", 16", 18" |

*Actual driving range will vary based on driving and battery conditions.
†Refer to Owner Manual for safe driving practices.